

DotDiver

GAI
Generic Assays





Compact, integrated system for **automated Immunodot Assays**

Unique Features

Flexible

Simultaneous performance of up to 24 different assays in an unprecedented short time.

Easy

Simply load the Immunodot Strips, the reagent cartridges and the patient samples, press on and follow the test processing data through a user friendly touch screen.

Automatic

The tests are carried out fully automatically according to established protocols. This will not bind lab capacities.

Fast

The DotDiver is much faster than comparable manual test execution.

Save


Through a bar-code association between the strips and the reagent cartridges.


Reliable

No pumps - no dead volume - no liquid handling
The standardized, automated test process ensures repeatable results.

Space saving

As compact as a microwave oven.





Comfortable position for sample application



- ○ ● **flexible**
- ○ ● **easy**
- ○ ● **automatic**

Consumables

Ready to use cartridges



Strips



Test properties

By repeatedly diving there is an optimal reaction environment, and the incubation time decreases. Specified protocols control the entire testing process, identical for all strips.

Accessories

Scanner for strips

With the optional scanner, the results can be automatically read from the strips and digitalized.

Analysis software: Dr.Dot

This special designed software enables a semiquantitative or quantitative (Quantrix assay) evaluation of the individual test results.

