

where treatment begins...

Pricing Valid to 31/03/2022

Date:	<input type="text"/>	Tel:	<input type="text"/>
Name:	<input type="text"/>	Fax:	<input type="text"/>
Invoice Address:	<input type="text"/>	E-mail:	<input type="text"/>
P O Number:	<input type="text"/>	Signature:	<input type="text"/>

"DuoTip II" is the most cost effective SPT system available. Uses just 5 microlitres of reagent per SPT. The Dipwell Tray comprises 40 Wells, each holding 0.25ml of each SPT allergen. Each individual DuoTip II lancet sits in the well of reagent and when lifted out the two tips of each lancet hold a meniscus of just 5 microlitres of the glycerinated SPT reagent. The circular scratch of the lancet introduces the very small volume of reagent into the epidermis. The twist-technique needs to be demonstrated and practised to get the optimal pressure and action, but is then easy for the Nurse or Doctor to use in clinical routine. Clinical documentation from formal evaluations is available on request, and on the Lincoln website at www.lincolndiagnostics.com

"UniTest PC" is the most recent innovative new lancet system from the USA's undisputed market leader. Uses the same Dipwell Tray and Wells as DuoTip, and the lancets cost the same. However the reagent drop size is approx. 10 microlitres, so reagent consumption is slightly greater than with DuoTip, but still far less than using a standard metal lancet with a normal drop size of 40-50 microlitres. The UniTest PC Lancet has the identical tip to the MultiTest PC Device and so gives exactly comparable clinical results. Clinical documentation from formal evaluations is available on request, and on the Lincoln website at www.lincolndiagnostics.com

"MultiTest PC" is the 3rd generation of the original MultiTest from the USA's market leader. Lincoln Diagnostics innovates, whilst other manufacturers merely try to emulate. PC stands for Pain Control. The tip of each of the 8 arms of the MultiTest PC device comprises a central pillar surrounded by 6 pointed tines. There are 3 devices sitting in the integral wells of the Tray, each well containing a known SPT allergen solution. Lifting the device from the Tray means the tines on each arm hold approx. 10 microlitres of SPT reagent. The device is pressed firmly into the patients forearm and gently rocked, then removed, and discarded. The central pillars prevent too deep penetration of the reagent-laden tines. The patient will notice that the degree of discomfort of the MultiTest PC device is considerably less than 8 individual SPTs with a metal lancet. The sole limitation of the MultiTest PC is the fact that fixed SPT panels of 8 allergens must be used. However, if any additional individual allergens require testing then use the UniTest PC which has identical lancets and so gives identical clinical responses. With MultiTest PC then 3 devices x 8 lancets = 24 allergens can be tested in under 30 seconds. Used together UniTest PC and MultiTest PC are the ultimate SPT system. Clinical documentation from formal evaluations is available on request, and on the Lincoln website at www.lincolndiagnostics.com

Other Lincoln Products

There are various other products manufactured by Lincoln, which we can obtain and make available for you. However the three products listed above, with their various components and accessories, are we know from international experience the most suitable and beneficial for UK consumers. However if you wish to receive any other Lincoln products than those stated here (as seen for example in the Lincoln website at www.lincolndiagnostics.com) then please contact us for product info and cost quotes.

Delivery

Delivery of all stated products is from local stock at Bio-Diagnostics. The delivery will then take 1-2 days by overnight courier. However if for any reason a product is not available from local stock then it will be ordered and imported from Lincoln Diagnostics of USA taking maximum 3-4 weeks. Please specify the recipient carefully, including the delivery address.

Prices

- ✓ Delivery cost is **£14.95** (ex-VAT) **£39.50** Ireland
- ✓ All prices exclude UK VAT at **20%**.
- ✓ Prices valid to **31/03/2022**
- ✓ Note all prices are eligible for negotiation, based primarily on volume commitment.

Payment Terms

- ✓ Payment needs to be made by Credit Card or Debit Card upon placement of this order with Bio-Diagnostics, or on account.
- ✓ Account Number: _____
- ✓ Contact Bio-Diagnostics to set up an account.

Delivery Address: _____

Lincoln Diagnostics Lancet Products

Product	SKU	Description	Price	Order
Dipwell Tray, DuoTip-Test II	DWTII-1	For 40 Lancets	£46.00	
Dipwell Tray Wells	DTW-40	40 Wells	£22.50	
DuoTip-Test II Lancets	LDT-II 400	Box of 20 x 20 = 400 Lancets	£85.00	
UniTest PC Lancets	UT-400	Box of 20 x 20 = 400 Lancets	£85.00	
MultiTest PC Tray	MHDWC-1	For 3 MultiTest Devices = 24 SPT	£28.50	
MultiTest PC Devices	MH-6 PC	Box of 6 MultiTest PC Devices	£18.00	
MultiTest PC Devices	MH-BLK PC	Case of 100 x 6 MultiTest PC Devices	£1550.00	

Send Order to Bio-Diagnostics:

Fax: 01684 592501

Tel: 01684 592262

E-mail: enquiries@bio-diagnostics.co.uk

Mail: Bio-Diagnostics Ltd, Rectory Rd, Upton Industrial Estate, Upton-upon-Severn, Worcs, WR8 0LX, UK.